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Bibliografia.

- Greenman, C., et al., Patterns of somatic mutation in human cancer genomes. *Nature*, 2007. 446(7132): p.153-8.
- Stephens, P.J., et al., The landscape of cancer genes and mutational processes in breast cancer. *Nature*, 2012. 486(7403): p. 400-4.
- Diaz, L.A., Jr. and A. Bardelli, Liquid biopsies: genotyping circulating tumor DNA. *J Clin Oncol*, 2014. 32(6): p. 579-86.
- Heitzer, E., P. Ulz, and J.B. Geigl, Circulating tumor DNA as a liquid biopsy for cancer. *Clin Chem*, 2015. 61(1): p. 112-23.
- Lebofsky, R., et al., Circulating tumor DNA as a non-invasive substitute to metastasis biopsy for tumor genotyping and personalized medicine in a prospective trial across all tumor types. *Mol Oncol*, 2015. 9(4): p. 783-90.
- Esposito, A., et al., Monitoring tumor-derived cell-free DNA in patients with solid tumors: clinical perspectives and research opportunities. *Cancer Treat Rev*, 2014. 40(5): p. 648-55.
- Newman, A.M., et al., An ultrasensitive method for quantitating circulating tumor DNA with broad patient coverage. *Nat Med*, 2014. 20(5): p. 548-54.
- Kidess, E., et al., Mutation profiling of tumor DNA from plasma and tumor tissue of colorectal cancer patients with a novel, high-sensitivity multiplexed mutation detection platform. *Oncotarget*, 2015. 6(4): p. 2549-61.
- Forshew, T., et al., Noninvasive identification and monitoring of cancer mutations by targeted deep sequencing of plasma DNA. *Sci Transl Med*, 2012. 4(136): p. 136ra68.
- Siegel, R.L., K.D. Miller, and A. Jemal, Cancer statistics, 2015. *CA Cancer J Clin*, 2015. 65(1): p. 5-29.
- Davies, H., et al., Mutation of the BRAF gene in human cancer. *Nature*, 2002. 417(6892): p. 949-54.
- Romero, A., et al., Identification of E545k mutation in plasma from a PIK3CA wild-type metastatic breast cancer patient by array-based digital polymerase chain reaction: Circulating-free DNA a powerful tool for biomarker testing in advance disease. *Transl Res*, 2015.
- Ascierto, P.A., et al., Phase II trial (BREAK-2) of the BRAF inhibitor dabrafenib (GSK2118436) in patients with metastatic melanoma. *J Clin Oncol*, 2013. 31(26): p. 3205-11.
- Rothschild, S.I., Targeted Therapies in Non-Small Cell Lung Cancer-Beyond EGFR and ALK. *Cancers (Basel)*, 2015. 7(2): p. 930-49.
- Heidary, M., et al., The dynamic range of circulating tumor DNA in metastatic breast cancer. *Breast Cancer Res*, 2014. 16(4): p. 421.
- Zill, O.A., et al., Cell-Free DNA Next-Generation Sequencing in Pancreatobiliary Carcinomas. *Cancer Discov*, 2015.
- Janku, F., et al., Actionable mutations in plasma cell-free DNA in patients with advanced cancers referred for experimental targeted therapies. *Oncotarget*, 2015. 6(14): p. 12809-21.
- Tabernero, J., et al., Analysis of circulating DNA and protein biomarkers to predict the clinical activity of regorafenib and assess prognosis in patients with metastatic colorectal cancer: a retrospective, exploratory analysis of the CORRECT trial. *Lancet Oncol*, 2015.
- Forbes, S.A., et al., COSMIC: exploring the world's knowledge of somatic mutations in human cancer. *Nucleic Acids Res*, 2015. 43(Database issue): p. D805-11.
- Shah, S.P., et al., Mutation of FOXL2 in granulosa-cell tumors of the ovary. *N Engl J Med*, 2009. 360(26): p.2719-29.
- Schirripa, M., et al., Role of NRAS mutations as prognostic and predictive markers in metastatic colorectal cancer. *Int J Cancer*, 2015. 136(1): p. 83-90.
- Janku, F., et al., PIK3CA mutations frequently coexist with RAS and BRAF mutations in patients with advanced cancers. *PLoS One*, 2011. 6(7): p. e22769.
- Kalfa, N., et al., Activating mutations of Gsalpha in kidney cancer. *J Urol*, 2006. 176(3): p. 891-5.
- Fecteau, R.E., et al., GNAS mutations identify a set of right-sided, RAS mutant, villous colon cancers. *PLoS One*, 2014. 9(1): p. e87966.
- Sparks, A.B., et al., Mutational analysis of the APC/beta-catenin/Tcf pathway in colorectal cancer. *Cancer Res*, 1998. 58(6): p. 1130-4.
- Bettegowda, C., et al., Detection of circulating tumor DNA in early- and late-stage human malignancies. *Sci Transl Med*, 2014. 6(224): p. 224ra24.
- Oshiro, C., et al., PIK3CA mutations in serum DNA are predictive of recurrence in primary breast cancer patients. *Breast Cancer Res Treat*, 2015. 150(2): p. 299-307.
- Beaver, J.A., et al., Detection of cancer DNA in plasma of patients with early-stage breast cancer. *Clin Cancer Res*, 2014. 20(10): p. 2643-50.
- Dawson, S.J., et al., Analysis of circulating tumor DNA to monitor metastatic breast cancer. *N Engl J Med*, 2013. 368(13): p. 1199-209.
- Ignatiadis, M. and S.J. Dawson, Circulating tumor cells and circulating tumor DNA for precision medicine: dream or reality? *Ann Oncol*, 2014. 25(12): p. 2304-13.
- Murtaza, M., et al., Non-invasive analysis of acquired resistance to cancer therapy by sequencing of plasma DNA. *Nature*, 2013. 497(7447): p. 108-12.
- Mohamed Suhaimi, N.A., et al., Non-invasive sensitive detection of KRAS and BRAF mutation in circulating tumor cells of colorectal cancer patients. *Mol Oncol*, 2015.

33. Sanmamed, M.F., et al., Quantitative cell-free circulating BRAFV600E mutation analysis by use of droplet digital PCR in the follow-up of patients with melanoma being treated with BRAF inhibitors. *Clin Chem*, 2015;61(1):p.297-304.
34. Siravegna, G., et al., Clonal evolution and resistance to EGFR blockade in the blood of colorectal cancer patients. *Nat Med*, 2015.
35. Roschewski, M., et al., Circulating tumour DNA and CT monitoring in patients with untreated diffuse large B-cell lymphoma: a correlative biomarker study. *Lancet Oncol*, 2015.
36. Heitzer E, Ulz P, Geigl JB. Circulating tumor DNA as a liquid biopsy for cancer. *Clinical chemistry*. 2015;61:112-23.
37. Lebofsky R, Decraene C, Bernard V, et al. Circulating tumor DNA as a non-invasive substitute to metastasis biopsy for tumor genotyping and personalized medicine in a prospective trial across all tumor types. *Molecular oncology*. 2015;9:783-90.
38. Esposito A, Bardelli A, Criscitiello C, et al. Monitoring tumor-derived cell-free DNA in patients with solid tumors: clinical perspectives and research opportunities. *Cancer treatment reviews*. 2014;40:648-55.
39. Diaz LA Jr, Bardelli A. Liquid biopsies: genotyping circulating tumor DNA. *Journal of clinical oncology : official journal of the American Society of Clinical Oncology*. 2014;32:579-86.
40. Perrone F, Lampis A, Bertan C, et al. Circulating free DNA in a screening program for early colorectal cancer detection. *Tumori*. 2014;100:115-21.
41. Amberger JS, Bocchini CA, Schiettecatte F, Scott AF, Hamosh A. OMIM.org: Online Mendelian Inheritance in Man (OMIM®), an online catalog of human genes and genetic disorders. *Nucleic acids research*. 2015;43:D789-98.
42. Forbes SA, Beare D, Gunasekaran P, et al. COSMIC: exploring the world's knowledge of somatic mutations in human cancer. *Nucleic acids research*. 2015;43:D805-11.
43. Ascierto PA, Minor D, Ribas A, et al. Phase II trial (BREAK-2) of the BRAF inhibitor dabrafenib (GSK2118436) in patients with metastatic melanoma. *Journal of clinical oncology : official journal of the American Society of Clinical Oncology*. 2013;31:3205-11.
44. Smalley KS, Xiao M, Villanueva J, et al. CRAF inhibition induces apoptosis in melanoma cells with non-V600E BRAF mutations. *Oncogene*. 2009;28:85-94.
45. Homet Moreno B, Ribas A. Anti-programmed cell death protein-1/ligand-1 therapy in different cancers. *British journal of cancer*. 2015;112:1421-7.
46. Morelli MP, Overman MJ, Dasari A, et al. Characterizing the patterns of clonal selection in circulating tumor DNA from patients with colorectal cancer refractory to anti-EGFR treatment. *Annals of oncology : official journal of the European Society for Medical Oncology/ESMO*. 2015;26:731-6.
47. Zill OA, Greene C, Sebisanoovic D, et al. Cell-Free DNA Next-Generation Sequencing in Pancreatobiliary Carcinomas. *Cancer discovery*. 2015;
48. Thress KS, Pawelczak CP, Felip E, et al. Acquired EGFR C797S mutation mediates resistance to AZD9291 in non-small cell lung cancer harboring EGFR T790M. *Nature medicine*. 2015;21:560-2.
49. Karachaliou N, Mayo-de Las Casas C, Queralt C, et al. Association of EGFR L858R Mutation in Circulating Free DNA With Survival in the EURTAC Trial. *JAMA oncology*. 2015;1:149-57.
50. Shah SP, Köbel M, Senz J, et al. Mutation of FOXL2 in granulosa-cell tumors of the ovary. *The New England journal of medicine*. 2009;360:2719-29.
51. De Stefano A, Carlomagno C. Beyond KRAS: Predictive factors of the efficacy of anti-EGFR monoclonal antibodies in the treatment of metastatic colorectal cancer. *World journal of gastroenterology : WJG*. 2014;20:9732-43.
52. Siravegna G, Mussolin B, Buscarino M, et al. Clonal evolution and resistance to EGFR blockade in the blood of colorectal cancer patients. *Nature medicine*. 2015;21:795-801.
53. Schirripa M, Cremolini C, Loupakis F, et al. Role of NRAS mutations as prognostic and predictive markers in metastatic colorectal cancer. *International journal of cancer. Journal international du cancer*. 2015;136:83-90.
54. Wong AL, Lim JS, Sinha A, et al. Tumour pharmacodynamics and circulating cell free DNA in patients with refractory colorectal carcinoma treated with regorafenib. *Journal of translational medicine*. 2015;13:57.
55. Rodon J, Braña I, Siu LL, et al. Phase I dose-escalation and -expansion study of buparlisib (BKM120), an oral pan-Class I PI3K inhibitor, in patients with advanced solid tumors. *Investigational new drugs*. 2014;32:670-81.
56. Forshew T, Murfaza M, Parkinson C, et al. Noninvasive identification and monitoring of cancer mutations by targeted deep sequencing of plasma DNA. *Science translational medicine*. 2012;4:136ra68.
57. Madic J, Kiialainen A, Bidard FC, et al. Circulating tumor DNA and circulating tumor cells in metastatic triple negative breast cancer patients. *International journal of cancer. Journal international du cancer*. 2015;136:2158-65.
58. Bettgowda C, Sausen M, Leary RJ, et al. Detection of circulating tumor DNA in early- and late-stage human malignancies. *Science translational medicine*. 2014;6:224ra24.
59. NCCN Clinical Practice Guidelines in Oncology (www.nccn.org/professionals/physician_gls/f_guidelines.asp#site) Accessed 15 June 2015.
60. Van Cutsem E, Cervantes A, Nordlinger B, Arnold D; ESMO Guidelines Working Group. Metastatic colorectal cancer: ESMO Clinical Practice Guidelines for diagnosis, treatment and follow-up. *Ann Oncol*. 2014;25 Suppl 3:iii1-9.
61. Higgins MJ, Baselga J (2011) Targeted therapies for breast cancer. *J Clin Invest* 121(10):3797-803.
62. Cancer Genome Atlas Research Network (2014) Comprehensive molecular profiling of lung adenocarcinoma. *Nature* 511(7511):543-50.
63. Swanton C, Futreal A, Eisen T (2006) Her2-targeted therapies in non-small cell lung cancer. *Clin Cancer Res* 12(14 Pt 2):4377s- 4383s.
64. Nakamura H, Kawasaki N, Taguchi M, et al. (2005) Association of HER-2 overexpression with prognosis in nonsmall cell lung carcinoma: a metaanalysis. *Cancer* 103(9):1865-73.
65. Tan D, Deeb G, Wang J, et al. (2003) HER-2/neu protein expression and gene alteration in stage I-IIIa non-small-cell lung cancer: a study of 140 cases using a combination of high throughput tissue microarray, immunohistochemistry, and fluorescent in situ hybridization. *Diagn Mol Pathol* 12(4):201-11.
66. Slamon DJ, Leyland-Jones B, Shak S, et al. (2001) Use of chemotherapy plus a monoclonal antibody against HER2 for metastatic breast cancer that overexpresses HER2. *N Engl J Med* 344(11):783-92.

67. Bang YJ, Van Cutsem E, Feyereislova A, et al. (2010) Trastuzumab in combination with chemotherapy versus chemotherapy alone for treatment of HER2-positive advanced gastric or gastro-oesophageal junction cancer (ToGA): a phase 3, open-label, randomised controlled trial. *Lancet* 376(9742):687-97.
68. Chumsri S, Weidler J, Ali S, et al. (2015) Prolonged Response to Trastuzumab in a Patient With HER2-Nonamplified Breast Cancer With Elevated HER2 Dimerization Harboring an ERBB2 S310F Mutation. *J Natl Compr Canc Netw* 13(9):1066-70.
69. Cappuzzo F, Bemis L, Varella-Garcia M (2006) HER2 mutation and response to trastuzumab therapy in non-small-cell lung cancer. *N Engl J Med* 354(24):2619-21.
70. Falchook GS, Janku F, Tsao AS, et al. (2013) Non-small-cell lung cancer with HER2 exon 20 mutation: regression with dual HER2 inhibition and anti-VEGF combination treatment. *J Thorac Oncol* 8(2):e19-20.
71. Mazières J, Peters S, Lepage B, et al. (2013) Lung cancer that harbors an HER2 mutation: epidemiologic characteristics and therapeutic perspectives. *J Clin Oncol* 31(16):1997-2003.
72. Baselga J, Cortés J, Kim SB, et al. (2012) Pertuzumab plus trastuzumab plus docetaxel for metastatic breast cancer. *N Engl J Med* 366(2):109-19.
73. Swain SM, Baselga J, Kim SB, et al. (2015) Pertuzumab, trastuzumab, and docetaxel in HER2-positive metastatic breast cancer. *N Engl J Med* 372(8):724-34.
74. Verma S, Miles D, Gianni L, et al. (2012) Trastuzumab emtansine for HER2-positive advanced breast cancer. *N Engl J Med* 367(19):1783-91.
75. Cameron D, Casey M, Oliva C, et al. (2010) Lapatinib plus capecitabine in women with HER-2-positive advanced breast cancer: final survival analysis of a phase III randomized trial. *Oncologist* 15(9):924-34.
76. Geyer CE, Forster J, Lindquist D, et al. (2006) Lapatinib plus capecitabine for HER2-positive advanced breast cancer. *N Engl J Med* 355(26):2733-43.
77. Serra V, Vivancos A, Puente XS, et al. (2013) Clinical response to a lapatinib-based therapy for a Li-Fraumeni syndrome patient with a novel HER2V659E mutation. *Cancer Discov* 3(11):1238-44.
78. Ali SM, Alpaugh RK, Downing SR, et al. (2014) Response of an ERBB2-Mutated Inflammatory Breast Carcinoma to Human Epidermal Growth Factor Receptor 2-Targeted Therapy. *J Clin Oncol ePub Feb 2014*.
79. Lin NU, Winer EP, Wheatley D, et al. (2012) A phase II study of afatinib (BIBW 2992), an irreversible ErbB family blocker, in patients with HER2-positive metastatic breast cancer progressing after trastuzumab. *Breast Cancer Res Treat* 133(3):1057-65.
80. Schwab CL, Bellone S, English DP, et al. (2014) Afatinib demonstrates remarkable activity against HER2-amplified uterine serous endometrial cancer in vitro and in vivo. *Br J Cancer* 111(9):1750-6.
81. De Grève J, Teugels E, Geers C, et al. (2012) Clinical activity of afatinib (BIBW 2992) in patients with lung adenocarcinoma with mutations in the kinase domain of HER2/neu. *Lung Cancer* 76(1):123-7.
82. De Grève J, Moran T, Graas MP, et al. (2015) Phase II study of afatinib, an irreversible ErbB family blocker, in demographically and genotypically defined lung adenocarcinoma. *Lung Cancer* 88(1):63-9.
83. Gandhi L, Bahleda R, Tolaney SM, et al. (2014) Phase I study of neratinib in combination with temsirolimus in patients with human epidermal growth factor receptor 2-dependent and other solid tumors. *J Clin Oncol* 32(2):68-75.
84. Ben-Baruch NE, Bose R, Kavuri SM, et al. (2015) HER2-Mutated Breast Cancer Responds to Treatment With Single-Agent Neratinib, a Second-Generation HER2/EGFR Tyrosine Kinase Inhibitor. *J Natl Compr Canc Netw* 13(9):1061-4.
85. Kris MG, Camidge DR, Giaccone G, et al. (2015) Targeting HER2 aberrations as actionable drivers in lung cancers: phase II trial of the pan-HER tyrosine kinase inhibitor dacomitinib in patients with HER2-mutant or amplified tumors. *Ann Oncol ePub Apr 2015*.
86. Takada M, Higuchi T, Tozuka K, et al. (2013) Alterations of the genes involved in the PI3K and estrogen-receptor pathways influence outcome in human epidermal growth factor receptor 2-positive and hormone receptor-positive breast cancer patients treated with trastuzumab-containing neoadjuvant chemotherapy. *BMC Cancer* 13:241.
87. Jensen JD, Knoop A, Laenkholm AV, et al. (2012) PIK3CA mutations, PTEN, and pHER2 expression and impact on outcome in HER2-positive early-stage breast cancer patients treated with adjuvant chemotherapy and trastuzumab. *Ann Oncol* 23(8):2034-42.
88. Berns K, Horlings HM, Hennessy BT, et al. (2007) A functional genetic approach identifies the PI3K pathway as a major determinant of trastuzumab resistance in breast cancer. *Cancer Cell* 12(4):395-402.
89. Dave B, Migliaccio I, Gutierrez MC, et al. (2011) Loss of phosphatase and tensin homolog or phosphoinositol-3 kinase activation and response to trastuzumab or lapatinib in human epidermal growth factor receptor 2-overexpressing locally advanced breast cancers. *J Clin Oncol* 29(2):166-73.
90. Loibl S, von Minckwitz G, Schneeweiss A, et al. (2014) PIK3CA Mutations Are Associated With Lower Rates of Pathologic Complete Response to Anti-Human Epidermal Growth Factor Receptor 2 (HER2) Therapy in Primary HER2-Overexpressing Breast Cancer. *J Clin Oncol ePub Sep 2014*.
91. Barbareschi M, Cuorvo LV, Girlando S, et al. (2012) PIK3CA mutations and/or PTEN loss in Her2-positive breast carcinomas treated with trastuzumab are not related to resistance to anti-Her2 therapy. *Virchows Arch* 461(2):129-39.
92. Guarnieri V, Generali DG, Frassoldati A, et al. (2014) Double-blind, placebo-controlled, multicenter, randomized, phase IIb neoadjuvant study of letrozole-lapatinib in postmenopausal hormone receptor-positive, human epidermal growth factor receptor 2-negative, operable breast cancer. *J Clin Oncol* 32(10):1050-7.
93. Jones KL, Buzdar AU (2009) Evolving novel anti-HER2 strategies. *Lancet Oncol* 10(12):1179-87.
94. Zagouri F, Sergentanis TN, Chrysikos D, et al. (2013) Hsp90 inhibitors in breast cancer: a systematic review. *Breast* 22(5):569-78.
95. Aarthy, R., Mani, S., Velusami, S. et al. *Mol Diagn Ther* (2015) 19: 339. doi:10.1007/s40291-015-0167-y

97. Chaudhuri AA, Binkley MS, Osmundson EC, Alizadeh AA, Diehn M. Predicting radiotherapy responses and treatment outcomes through analysis of circulating tumor DNA. *Seminars in radiation oncology*. 2015;25(4):305-312. doi:10.1016/j.semradonc.2015.05.001.
98. Michail Ignatiadis, Mark Lee and Stefanie S. Jeffrey, *Clin Cancer Res* November 1 2015 (21) (21) 4786-4800; DOI: 10.1158/1078-0432.CCR-14-1190
99. Jiri Polivka Jr, Martin Pesta, and Filip Janku, *Expert Review Of Molecular Diagnostics* Vol. 15 , Iss. 12,2015
100. Diehl F, Schmidt K, Choti MA, et al. Circulating mutant DNA to assess tumor dynamics. *Nature medicine*. 2008;14(9):985-990. doi:10.1038/nm.1789.
101. Diehl F, Li M, Dressman D, et al. Detection and quantification of mutations in the plasma of patients with colorectal tumors. *Proceedings of the National Academy of Sciences of the United States of America*. 2005;102(45):16368-16373. doi:10.1073/pnas.0507904102.
102. Fan HC, Blumenfeld YJ, Chitkara U, Hudgins L, Quake SR. Noninvasive diagnosis of fetal aneuploidy by shotgun sequencing DNA from maternal blood. *Proceedings of the National Academy of Sciences of the United States of America*. 2008;105(42):16266-16271. doi:10.1073/pnas.0808319105.
103. Sabine Jahr, Hannes Hentze, Sabine Englisch, Dieter Hardt, Frank O. Fackelmayer, Rolf-Dieter Heschand Rolf Knippers, *Cancer Res* February 2 2001 (61) (4) 1659-1665
104. Mouliere F, Robert B, Arnau Peyrotte E, et al. High Fragmentation Characterizes Tumour-Derived Circulating DNA. Lee T, ed. *PLoS ONE*. 2011;6(9):e23418. doi:10.1371/journal.pone.0023418.
105. Holdhoff et al., *JNCI J Natl Cancer Inst* (2009) 101(18): 1284-1285. doi: 10.1093/jnci/djp240
106. Newman AM, Bratman SV, To J, et al. An ultrasensitive method for quantitating circulating tumor DNA with broad patient coverage. *Nature medicine*. 2014;20(5):548-554. doi:10.1038/nm.3519.
107. Cécile Jovelet et al., *Clin Cancer Res* June 15 2016 (22) (12) 2960-2968
108. Gerlinger M, Rowan AJ, Horswell S, et al. Intratumor Heterogeneity and Branched Evolution Revealed by Multiregion Sequencing. *The New England journal of medicine*. 2012;366(10):883-892. doi:10.1056/NEJMoa1113205.
109. Hiley C, de Bruin EC, McGranahan N, Swanton C. Deciphering intratumor heterogeneity and temporal acquisition of driver events to refine precision medicine. *Genome Biology*. 2014;15(8):453. doi:10.1186/s13059-014-0453-8.
110. Ichihara E, Lovly CM. Shades of T790M – intratumor heterogeneity in EGFR mutant lung cancer. *Cancer discovery*. 2015;5(7):694-696. doi:10.1158/2159-8290.CD-15-0616.
111. Nik-Zainal S, Van Loo P, Wedge DC, et al. The Life History of 21 Breast Cancers. *Cell*. 2012;149(5):994-1007. doi:10.1016/j.cell.2012.04.023.
112. Piotrowska Z, Niederst MJ, Karlovich CA, et al. Heterogeneity Underlies the Emergence of EGFR T790 Wild-Type Clones Following Treatment of T790M-Positive Cancers with a Third Generation EGFR Inhibitor. *Cancer discovery*. 2015;5(7):713-722. doi:10.1158/2159-8290.CD-15-0399.
113. Wang Y, Waters J, Leung ML, et al. Clonal Evolution in Breast Cancer Revealed by Single Nucleus Genome Sequencing. *Nature*. 2014;512(7513):155-160. doi:10.1038/nature13600.
114. Patel eTsui, 2015, *Clinical Biochemistry*, Volume 48, Issue 15, October 2015, Pages 957–961